


2022
product catalog

GERMAN INDUSTRY
CE / ISO 9001



■ GYNECOLOGY

Developed and manufactured in Germany

 RG Medical USA, Inc. |  15065 Lebanon Road, Suite 201 Old Hickory, TN 37138 USA
 615-269-7256 ·  615-269-4605

www.rgendo.com

■ Medical monitors



- ▷ Screen size 23.8 "
- ▷ Resolution 1920 x 1080
- ▷ Brightness 600cd / m '
- ▷ Color depth
- ▷ 16.7 million colors
- ▷ Viewing angle 78/178
- ▷ Response time (GtoG) 14ms

M-238CL1

■ Cameras

1080p wireless camera



- ▷ Without USB
- ▷ Free software
- ▷ Images / Recording

EOC-23

■ Digital insufflator



5 AÑOS GARANTÍA

- ▷ 40 liter digital insufflator
- ▷ Touch screen
- ▷ Day and night mode

E-IN40

1080p USB Camera + 10W LED



- ▷ With USB
- ▷ Free software
- ▷ Images / Recording
- ▷ LED

EC-USB1
EC-USB1 L

■ Combo camera Full HD

Full HD + 300w camera. LED light source Records in Full HD 1080p 30 frames



- ▷ LED 300
- ▷ 1/3 "CMOS sensor
- ▷ 3.0 megapixels
- ▷ Color: 0.1Lux@F1.2
- ▷ 1920x1080 / 60Hz Full HD
- ▷ 3G SDI 1080p
- ▷ HDMI 1080p
- ▷ Composite video C-Video
- ▷ 5300 K

E-CHD11

* The set includes camera, universal cable and equipment.

■ Passive work item

Endoscopy Ø 2.9 mm



<input checked="" type="checkbox"/>	17/18 Fr.	U-1- 530-17
	1 outer tube	U-1- 534-17
	1 inner tube	U-1- 532-17
	1 shutter	U-1- 580-17

<input type="checkbox"/>	19/22 Fr.	U-1- 530-19
	1 outer tube	U-1- 534-19
	1 inner tube	U-1- 532-19
	1 shutter	U-1- 580-19

Endoscopy Ø 4.0 mm



<input checked="" type="checkbox"/>	22/24Fr.	U-1- 530-22
	1 outer tube	U-1- 534-22
	1 inner tube	U-1- 532-22
	1 shutter	U-1- 580-22

<input checked="" type="checkbox"/>	24/26 Fr.	U-1- 530-24
	1 outer tube	U-1- 534-24
	1 inner tube	U-1- 532-24
	1 shutter	U-1- 580-24

■ Electrodes

Cut loop, retrograde 30°



<input checked="" type="checkbox"/>	17 Ch./Fr.	U-1- 200-17
<input type="checkbox"/>	19 Ch./Fr.	U-1- 200-19
<input checked="" type="checkbox"/>	22 Ch./Fr.	U-1- 200-22
<input checked="" type="checkbox"/>	24 Ch./Fr.	U-1- 200-24

Straight cut bow



<input checked="" type="checkbox"/>	17 Ch./Fr.	U-1- 202-17
<input type="checkbox"/>	19 Ch./Fr.	U-1- 202-19
<input checked="" type="checkbox"/>	22 Ch./Fr.	U-1- 202-22
<input checked="" type="checkbox"/>	24 Ch./Fr.	U-1- 202-24

Retrograde cut band



<input checked="" type="checkbox"/>	17 Ch./Fr.	U-1- 204-17
<input type="checkbox"/>	19 Ch./Fr.	U-1- 204-19
<input checked="" type="checkbox"/>	22 Ch./Fr.	U-1- 204-22
<input checked="" type="checkbox"/>	24 Ch./Fr.	U-1- 204-24

Knife electrode



<input checked="" type="checkbox"/>	17 Ch./Fr.	U-1- 206-17
<input type="checkbox"/>	19 Ch./Fr.	U-1- 206-19
<input checked="" type="checkbox"/>	22 Ch./Fr.	U-1- 206-22
<input checked="" type="checkbox"/>	24 Ch./Fr.	U-1- 206-24

Ball electrode
3mm roller



<input checked="" type="checkbox"/>	17 Ch./Fr.	U-1- 208-17
<input type="checkbox"/>	19 Ch./Fr.	U-1- 208-19
<input checked="" type="checkbox"/>	22 Ch./Fr.	U-1- 208-22
<input checked="" type="checkbox"/>	24 Ch./Fr.	U-1- 208-24

Ball electrode
5mm roller




<input checked="" type="checkbox"/>	17 Ch./Fr.	U-1- 209-17
<input type="checkbox"/>	19 Ch./Fr.	U-1- 209-19
<input checked="" type="checkbox"/>	22 Ch./Fr.	U-1- 209-22
<input checked="" type="checkbox"/>	24 Ch./Fr.	U-1- 209-24

■ Autoclave, 134° C / 273° F



<u>Vision angle</u>	<u>Ø 2.0 mm</u>	<u>Ø 2.9 mm</u>	<u>Ø 4.0 mm</u>
0°	E020-240-00	E029-303-00	E040-303-00
12°	E020-240-12	E029-303-12	E040-303-12
30°	E020-240-30	E029-303-30	E040-303-30
70°	E020-240-70	E029-303-70	E040-303-70

■ Diagnostic shirt for hysteroscopy



<u>Scopes</u>	<u>Ext. sheath</u>	<u>Diag. Sheath</u>	<u>Obturador</u>
Ø 2.9 mm	30-022-05	30-022-06	30-022-00
Ø 4.0 mm	30-023-05	30-023-06	30-023-00

■ Camisa operatoria de Bettocchi "Style"



<u>scopio</u>	<u>Working channel Ø</u>	
Ø 2.9 mm	5 Charr./Fr	30-018-05
Ø 4.0 mm	7 Charr./Fr.	30-018-07

■ Forceps



No insulated
 Ø 7 Ch./Fr . 30-055-07
 Scissor Ø 5 Ch./Fr . 30-055-05

No insulated
 Ø 7 Ch./Fr . 30-050-07
 Biopsy Clamp Ø 5 Ch./Fr . 30-050-05

No insulated
 Ø 7 Ch./Fr . 30-051-07
 Alligator Clamp Ø 5 Ch./Fr . 30-051-05



Tischler 30-070-01
 Baby Tischler 30-070-02


Kevorkian 30-070-03

■ Electrodes




Loop electrode


5 mm x 5 mm M90-001




10 mm x 5 mm M90-002
 10 mm x 7 mm M90-003
 10 mm x 10 mm M90-004



15 mm x 5 mm M90-011
 15 mm x 7 mm M90-012
 15 mm x 8 mm M90-013
 15 mm x 15 mm M90-014



20 mm x 8 mm M90-021
 20 mm x 10 mm M90-022
 20 mm x 13 mm M90-023
 20 mm x 15 mm M90-024
 20 mm x 20 mm M90-025
 20 mm x 20 mm M90-026
 25 mm x 10 mm M90-027



Square Electrode

10 mm x 4 mm M90-031
 10 mm x 8 mm M90-032



Ball Electrode

3 mm M90-051
 5 mm M90-052



■ Light cable



Size

Length 180 cm E1-035-18

■ Video Colposcope



Video Colposcope

M99-001

■ Mini 100W LED light source



- 100W LED
- German design quality
- 20.000 + hours
- Auto-adjusting light intensity

**+ 60%
Light**

Specifications

Light Type	100 Watt LED
Ports	Single (Storz), other adapters optional.
Light Intensity Control	Auto Adjust
Color Temp.	6000-6500°K
LED Array Life	20.000 + Hours
Dimensions	9x12x7cm
Weight	280 grs
Warranty	2 Years

Limited Warranty: RG Medical instruments are warranted to be free from defects in materials and workmanship for the life of the instrument. Any instrument found to be defective will be repaired or replaced at RG Medical's discretion. This warranty is void if the instrument is not used for the surgical purpose for which it was designed or if it is serviced by someone other than RG Medical.

* Please note that the manufacturer only provides a stock image for this product group. The appearance of the product may be different. Specific product details are listed in the description.

■ Portable 10W LED light



Features

Compact and portable, the built-in lithium battery allows the mini portable cold light source to work for more than 2 hours after recharging. We design and manufacture it for use with a rigid or flexible endoscope.

A particularly well-conceived safety feature is the automatic switch from the main light source to the spare light source in the event of a light source failure.

Therefore, the light source is always ready for emergency use and contributes to patient safety.



Includes:

- LED device
- 2 rechargeable batteries
- Charger
- Optional small travel case



10W LED head-light **E-HL10P**



GERMAN INDUSTRY
CE / ISO 9001