



>> Mini 100W Led light source

■ MINI 100W LED LIGHT SOURCE



- 100W LED
- German design quality
- 20.000 + hours
- Auto-adjusting light intensity

**+ 60%
Light**

Specifications

Light Type	100 Watt LED
Ports	Single (Storz), other adapters optional.
Light Intensity Control	Auto Adjust
Color Temp.	6000-6500°K
LED Array Life	20.000 + Hours
Dimensions	9x12x7cm
Weight	280 grs
Warranty	2 Years

Limited Warranty: RG Medical instruments are warranted to be free from defects in materials and workmanship for the life of the instrument. Any instrument found to be defective will be repaired or replaced at RG Medical's discretion. This warranty is void if the instrument is not used for the surgical purpose for which it was designed or if it is serviced by someone other than RG Medical.

* Please note that the manufacturer only provides a stock image for this product group. The appearance of the product may be different. Specific product details are listed in the description.

 RG Medical USA, Inc. |  15065 Lebanon Road, Suite 201 Old Hickory, TN 37138 USA

 615-269-7256 ·  615-269-4605

www.rgendo.com



>> Portable 10w Led light

■ PORTABLE 10W LED LIGHT



Compact and portable, the built-in lithium battery allows the mini portable cold light source to work for more than 2 hours after recharging.

We design and manufacture it for use with a rigid or flexible endoscope.

A particularly well-conceived safety feature is the automatic switch from the main light source to the spare light source in the event of a light source failure.

Therefore, the light source is always ready for emergency use and contributes to patient safety.



Way to use it in the optician.

Includes:

- ▶ LED device
- ▶ 2 rechargeable attitudes
- ▶ Charger
- ▶ Optional small travel case

10W LED head-light **E-HL10P**



 RG Medical USA, Inc. |  15065 Lebanon Road, Suite 201 Old Hickory, TN 37138 USA

 615-269-7256 ·  615-269-4605

www.rgendo.com

■ HEADLIGHT SPECIAL

10W ULTRA-BRIGHT & LITE Weight Headlight



*also with zoom option

10W ULTRA-BRIGHT & LITE Weight Headlight



LED 5W HEADLIGHT Wireless



Characteristics

Applications:
Surgical, for dentistry,
veterinary

Lighting type:
LED

Other features:
With rechargeable battery, for binocular
loupe, with portable light source

Headlights include:



- Portable travel case
- 1 Charger
- 2 rechargeable lithium-ion batteries