



>> Mini 100W Led light source

■ MINI 100W LED LIGHT SOURCE



■ **100W LED**

■ **GERMAN DESIGN QUALITY**

■ **20.000 + HOURS**

■ **AUTO-ADJUSTING LIGHT INTENSITY**

**+ 60%
Light**

Specifications

Light Type	100 Watt LED
Ports	Single (Storz), other adapters optional.
Light Intensity Control	Auto Adjust
Color Temp.	6000-6500°K
LED Array Life	20.000 + Hours
Dimensions	9x12x7cm
Weight	280 grs
Warranty	2 Years

Limited Warranty: RG Medical instruments are warranted to be free from defects in materials and workmanship for the life of the instrument. Any instrument found to be defective will be repaired or replaced at RG Medical's discretion. This warranty is void if the instrument is not used for the surgical purpose for which it was designed or if it is serviced by someone other than RG Medical.

* Please note that the manufacturer only provides a stock image for this product group. The appearance of the product may be different. Specific product details are listed in the description.

 RG Medical USA, Inc. |  15065 Lebanon Road, Suite 201 Old Hickory, TN 37138 USA

 615-269-7256 ·  615-269-4605

www.rgendo.com