

URETHERO-RENSCOPE FLEXIBLE



- No Cross Infection**
- No Maintenance**
- No Sterilization**
- No optical fiber**
- Handheld, portable, easy use**



Greater Tip Range of Motion



Passive Bending



Two-way Deflection




1:1 Torque Ratio

SPECIFICATION

Optical System	
Field of view	90°
Direction of vision	0° (forward)
Depth of view field	3-50 mm
Image Resolution	1080 p (400x400 native pixels)
Led	
Angulation range	Up 270° Down 270°
Outer diameter insertion tube	7.5Fr / 9.5 Fr
Working Length	650 mm
Channel	
Working Channel	3.6 Fr
Irrigation	30 mL/min

*For medical professionals only

270° ANGULATION



HIGH QUALITY OF IMAGE



Image Resolution
1080 p (400x400 native pixels)