


AAO SPECIAL 2024

OTOSCOPE 3/4mm



Non-autoclavable From

Otoscope 3mm / 110mm


0°	E030-110-00
30°	E030-110-30

Autoclavable From

Otoscope 4mm / 110mm

0°	E040-110-00
30°	E040-110-30

SINUSCOPE 2.7/4mm



Non-autoclavable From

Sinuscope 3mm / 175mm

0°	E030-175-00
30°	E030-175-30
70°	E030-175-70

Autoclavable From

Sinuscope 4mm / 175mm

0°	E040-175-00
30°	E040-175-30
45°	E040-175-45
70°	E040-175-70

LARYNGOSCOPE 6/8mm



From

Laryngoscope 6mm

70°	E060-185-70
90°	E060-185-90

Laryngoscope 8mm

70°	E080-185-70
90°	E080-185-90

REPAIRS

Rigid Endoscope Repairs

- // Level I Interior cleaning, polishing fiber/objective, rebond adhesive seals
- // Level II Level I plus replacement of 2 rod or achromat lens

Flexible Endoscope Repairs

- // Level I Bending rubber replacement
- // Level II Level I plus angulation adjustment

FREE REPAIR QUOTE

1 year warranty on all repairs