

2024

Produktkatalog



# ■ LARYNGOSKOPIE

*Entwickelt und hergestellt in Deutschland*




**03** Handles



**04** Round cup forceps with Handle L-001-01



**05** Round cup forceps with Handle AA



**06** Oval cup forceps with Handle L-001-01

- Oval cup forceps with Handle AA



**07** Suction Tube mountable on the handle

- Triangular with Handle L-001-01



**08** Triangular with Handle AA

**09** Grasping forceps with Handle L-001-01

**10** Grasping forceps with Handle AA



- Suction tube mountable on the handle

**11** Scissors with Handle L-001-01

**12** Scissors with Handle AA



- 13** Round cup forceps Kleinsasser with Suction channel
- 14** Vocal cord spreader, Needle holder, Shaft curved forceps
- 15** Grasping forceps
- 16** Cutting forceps, complete matt finish with suction channel
  - Grasping forceps fine serrated, complete matt finish
- 17** Grasping forceps with Laser-Suction channel
  - Rotatable - Laryngo instruments
- 18** Monopolar HF - instruments
  - Flexible grasping ans forceps dor laryngo-fiberscope
  - Grasping and biopsy forceps for Slendoscope

2024

Product catalog

[www.rgmed.com](http://www.rgmed.com)



**Art. N°**  
L-001-01



**Art. N°**  
L-001-02



**Art. N°**  
L-001-03

Round cup forceps with Handle L-001-01



**Shape with conical Tube  
Shaft 0 3,0 -2.0mm**

Shape with conical Tube Shaft 0 3,0 -2.0mm	Lenght	Diameter	Art. N°
Round cup, jaw straight	230 mm	1,0 mm	L-001-04S
Round cup, jaw upwards	230 mm	1,0 mm	L-001-04U
Round cup, jaw right angled	230 mm	1,0 mm	L-001-04R
Round cup, jaw left angled	230 mm	1,0 mm	L-001-04L
Round cup, jaw straight	180 mm	1,0 mm	L-002-04S
Round cup, jaw upwards	180 mm	1,0 mm	L-002-04U
Round cup, jaw right angled	180 mm	1,0 mm	L-002-04R
Round cup, jaw left angled	180 mm	1,0 mm	L-002-04L
Round cup, jaw straight	250 mm	2,0 mm	L-003-05S
Round cup, jaw upwards	250 mm	2,0 mm	L-003-05U
Round cup, jaw right angled	250 mm	2,0 mm	L-003-05R
Round cup, jaw left angled	250 mm	2,0 mm	L-003-05L
Round cup, jaw straight	230 mm	3,0 mm	L-004-06S
Round cup, jaw upwards	230 mm	3,0 mm	L-004-06U
Round cup, jaw right angled	230 mm	3,0 mm	L-004-06R
Round cup, jaw left angled	230 mm	3,0 mm	L-004-06L
Round cup, jaw straight	230 mm	4,0 mm	L-005-07S
Round cup, jaw upwards	230 mm	4,0 mm	L-005-07U
Round cup, jaw right angled	230 mm	4,0 mm	L-005-07R
Round cup, jaw left angled	230 mm	4,0 mm	L-005-07L
Round cup, jaw straight	230 mm	5,0 mm	L-006-08S

### Round cup forceps with Handle AA



Shape with conical Tube Shaft 0 3,0 -2.0mm	Lenght	Diameter	Art. N°
Round cup, jaw straight	230 mm	1,0 mm	AA10-1S
Round cup, jaw upwards	230 mm	1,0 mm	AA10-1U
Round cup, jaw right angled	230 mm	1,0 mm	AA10-1R
Round cup, jaw left angled	230 mm	1,0 mm	AA10-1L
Round cup, jaw straight	180 mm	1,0 mm	AA15-1S
Round cup, jaw upwards	180 mm	1,0 mm	AA15-1U
Round cup, jaw right angled	180 mm	1,0 mm	AA15-1R
Round cup, jaw left angled	180 mm	1,0 mm	AA15-1L
Round cup, jaw straight	250 mm	2,0 mm	AA10-2S
Round cup, jaw upwards	250 mm	2,0 mm	AA10-2U
Round cup, jaw right angled	250 mm	2,0 mm	AA10-2R
Round cup, jaw left angled	250 mm	2,0 mm	AA10-2L
Round cup, jaw straight	230 mm	3,0 mm	AA10-3S
Round cup, jaw upwards	230 mm	3,0 mm	AA10-3U
Round cup, jaw right angled	230 mm	3,0 mm	AA10-3R
Round cup, jaw left angled	230 mm	3,0 mm	AA10-3L
Round cup, jaw straight	230 mm	4,0 mm	AA10-4S
Round cup, jaw upwards	230 mm	4,0 mm	AA10-4U
Round cup, jaw right angled	230 mm	4,0 mm	AA10-4R
Round cup, jaw left angled	230 mm	4,0 mm	AA10-4L
Round cup, jaw straight	230 mm	5,0 mm	AA10-5S

With conical shaft and irrigation port

Oval cup forceps with Handle L-001-01



Shape with conical Tube Shaft Ø 3,0 -2.0mm	Lenght	Diameter	Art. N°
Oval cup, jaw straight	230 mm	2,0 mm	LAR18-2S
Oval cup, jaw straight	230 mm	3,0 mm	LAR18-3S
Oval cup, jaw straight	230 mm	4,0 mm	LAR18-4S
Oval cup, jaw straight angeld	230 mm	2,0 mm	LAR18-2U
Oval cup, jaw straight angeld	230 mm	3,0 mm	LAR18-3U
Oval cup, jaw straight angeld	230 mm	4,0 mm	LAR18-4U

Oval cup forceps with Handle AA



Shape with conical Tube Shaft Ø 3,0 -2.0mm	Lenght	Diameter	Art. N°
Oval cup, jaw straight	230 mm	2,0 mm	AA18-2S
Oval cup, jaw straight	230 mm	3,0 mm	AA18-3S
Oval cup, jaw straight	230 mm	4,0 mm	AA18-4S
Oval cup, jaw straight angeld	230 mm	2,0 mm	AA18-2U
Oval cup, jaw straight angeld	230 mm	3,0 mm	AA18-3U
Oval cup, jaw straight angeld	230 mm	4,0 mm	AA18-4U


Suction Tube mountable on the handle



The suction tube can be ordered separately

Art. N°  
L-001-04

Triangular with Handle L-001-01



**Shape with conical Tube  
Shaft Ø 3,0 -2.0mm**

**Lenght**

**Diameter**

**Art. N°**

Triangular, jaw straight	230 mm	2,0 mm	LART10-2S
Triangular, jaw right angled	230 mm	2,0 mm	LART10-2R
Triangular, jaw left angled	230 mm	2,0 mm	LART10-2L
Triangular, jaw upwards	230 mm	2,0 mm	LART10-2U
Triangular, jaw straight	230 mm	3,0 mm	LART10-3S
Triangular, jaw right angled	230 mm	3,0 mm	LART10-3R
Triangular, jaw left angled	230 mm	3,0 mm	LART10-3L
Triangular, jaw upwards	230 mm	3,0 mm	LART10-3U



## Triangular with Handle AA



**Shape with conical Tube  
Shaft 0 3,0 -2.0mm**

**Lenght**


**Diameter**

**Art. N°**

Triangular, jaw straight	230 mm	2,0 mm	AA110-2S
Triangular, jaw right andled	230 mm	2,0 mm	AA110-2R
Triangular, jaw left angled	230 mm	2,0 mm	AA110-2L
Triangular, jaw upwards	230 mm	2,0 mm	AA110-2U
Triangular, jaw straight	230 mm	3,0 mm	AA110-3S
Triangular, jaw right angled	230 mm	3,0 mm	AA110-3R
Triangular, jaw left angled	230 mm	3,0 mm	AA110-3L
Triangular, jaw upwards	230 mm	3,0 mm	AA110-3U

Other Tube diameter upon request available

Grasping forceps with Handle L-001-01



Shape with conical Tube Shaft 0 3,0 -2.0mm	Lenght	Diameter	Art. N°
Straight, fine serrated	230 mm	1,0 mm	LAR20-1S
Left curved, fine serrated	230 mm	1,0 mm	LAR20-1L
Right curved, fine serrated	230 mm	1,0 mm	LAR20-1R
Straight, fine serrated	180 mm	1,0 mm	LAR25-1S
Left curved, fine serrated	180 mm	1,0 mm	LAR25-1L
Right curved, fine serrated	180 mm	1,0 mm	LAR25-1R
Straight jaw, fine serrated	230 mm	2,0 mm	LAR20-2S
Upwards curved jaw , fine serrated	230 mm	2,0 mm	LAR20-2U
Left curved, fine serrated	230 mm	2,0 mm	LAR20-2L
Right curved, fine serrated	230 mm	2,0 mm	LAR20-2R
Straight, fine serrated	180 mm	2,0 mm	LAR25-2S
Left curved, fine serrated	180 mm	2,0 mm	LAR25-2L
Right curved, fine serrated	180 mm	2,0 mm	LAR25-2R
Straight jaw, fine serrated	230 mm	3,0 mm	LAR20-3S
Upwards curved jaw , fine serrated	230 mm	3,0 mm	LAR20-3U
Left curved, fine serrated	230 mm	3,0 mm	LAR20-3L
Right curved, fine serrated	230 mm	3,0 mm	LAR20-3R
Straight, alligator teeth	230 mm	2,0 mm	LAR23-2S
Straight, alligator teeth	230 mm	3,0 mm	LAR23-3S

Other Tube diameter upon request available

## Grasping forceps with Handle AA



### Shape with conical Tube Shaft 0 3,0 -2.0mm

Shape with conical Tube Shaft 0 3,0 -2.0mm	Lenght	Diameter	Art. N°
Straight, fine serrated	230 mm	1,0 mm	AA20-1S
Left curved, fine serrated	230 mm	1,0 mm	AA20-1L
Right curved, fine serrated	230 mm	1,0 mm	AA20-1R
Straight, fine serrated	180 mm	1,0 mm	AA25-1S
Left curved, fine serrated	180 mm	1,0 mm	AA25-1L
Right curved, fine serrated	180 mm	1,0 mm	AA25-1R
Straight jaw, fine serrated	230 mm	2,0 mm	AA20-2S
Upwards curved jaw , fine serrated	230 mm	2,0 mm	AA20-2U
Left curved, fine serrated	230 mm	2,0 mm	AA20-2L
Right curved, fine serrated	230 mm	2,0 mm	AA20-2R
Straight, fine serrated	180 mm	2,0 mm	AA25-2S
Left curved, fine serrated	180 mm	2,0 mm	AA25-2L
Right curved, fine serrated	180 mm	2,0 mm	AA25-2R
Straight jaw, fine serrated	230 mm	3,0 mm	AA20-3S
Upwards curved jaw , fine serrated	230 mm	3,0 mm	AA20-3U
Left curved, fine serrated	230 mm	3,0 mm	AA20-3L
Right curved, fine serrated	230 mm	3,0 mm	AA20-3R
Straight, alligator teeth	230 mm	2,0 mm	AA23-2S
Straight, alligator teeth	230 mm	3,0 mm	AA23-3S


## Suction Tube mountable on the handle



The suction tube can be ordered separately

Art. N°  
L-001-05

Scissors with Handle L-001-01



**Shape with conical Tube  
Shaft Ø 3,0 -2.0mm**

Shape	Length	Diameter	Art. N°
Straight	230 mm	1,0 mm	LAR30-1S
Left curved	230 mm	1,0 mm	LAR30-1L
Right curved	230 mm	1,0 mm	LAR30-1R
Straight,	180 mm	1,0 mm	LAR35-1S
Left curved	180 mm	1,0 mm	LAR35-1L
Right curved	180 mm	1,0 mm	LAR35-1R
Straight	250 mm	2,0 mm	LAR30-2S
Left curved	250 mm	2,0 mm	LAR20-2L
Right curved	250 mm	2,0 mm	LAR20-2R
Straight	230 mm	3,0 mm	LAR30-3S
Left curved	230 mm	3,0 mm	LAR30-3L
Right curved	230 mm	3,0 mm	LAR30-3R



Shape	Length	Diameter	Art. N°
45° angled curved	230 mm	3,0 mm	LAR45-3U

Scissors with Handle AA



**Shape with conical Tube Shaft 0 3,0 -2.0mm**

	<b>Lenght</b>	<b>Diameter</b>	<b>Art. N°</b>
Straight	230 mm	1,0 mm	AA30-1S
Left curved	230 mm	1,0 mm	AA30-1L
Right curved	230 mm	1,0 mm	AA30-1R
Straight,	180 mm	1,0 mm	AA35-1S
Left curved	180 mm	1,0 mm	AA35-1L
Right curved	180 mm	1,0 mm	AA35-1R
Straight	250 mm	2,0 mm	AA30-2S
Left curved	250 mm	2,0 mm	AA20-2L
Right curved	250 mm	2,0 mm	AA20-2R
Straight	230 mm	3,0 mm	AA30-3S
Left curved	230 mm	3,0 mm	AA30-3L
Right curved	230 mm	3,0 mm	AA30-3R



<b>Shape</b>	<b>Lenght</b>	<b>Diameter</b>	<b>Art. N°</b>
45° angled curved	230 mm	3,0 mm	A45-3U


Round cup forceps Kleinsasser with Suction channel



**Shape with conical Tube  
Shaft 0,3,0 -2.0mm**

	<b>Length</b>	<b>Diameter</b>	<b>Art. N°</b>
Round cup, Straight	230 mm	2,0 mm	LAR40-2S
Round cup, Upwards	230 mm	2,0 mm	LAR40-2U
Round cup, Left	230 mm	2,0 mm	LAR40-2L
Round cup, Right	230 mm	2,0 mm	LAR40-2R
Round cup, Straight	230 mm	4,0 mm	LAR40-4S
Round cup, Upwards	230 mm	4,0 mm	LAR40-4U
Round cup, Left	230 mm	4,0 mm	LAR40-4L
Round cup, Right	230 mm	4,0 mm	LAR40-4R

Vocal cord spreader, Needle holder, Shaft curved forceps



Art. N°  
LAR 95-20

Shape with conical Tube Shaft 0 3,0 -2.0mm	Lenght	Shaft Diameter	Art. N°
Vocal cord spreader diameter 4,5	250 mm	2,0 mm	LAR70-20
Larynx-, Alligator jaw shaft 80° downwards	260 mm	2,0 mm	LAR90-20
Larynx-forceps, round jaw 45° shaft 80° downwards	260 mm	2,0 mm	LAR96-20
Larynx-, Alligator jaw shaft 80° downwards	260 mm	2,0 mm	LAR95-20
Needle holder forceps with ratchet handle	250 mm	2,0 mm	LAR100-20



Art. N°  
LAR70-20

Grasping forceps



Art. N°  
LAR50-1



Art. N°  
LAR50-2



Art. N°  
LAR50-3



Art. N°  
LAR50-4



Art. N°  
LAR50-5

**Shape with conical Tube  
Shaft 0 3,0 -2.0mm**

Shape with conical Tube Shaft 0 3,0 -2.0mm	Lenght	Diameter	Art. N°
Fine serrated, small	230 mm	2,0 mm	LAR50-1
Alligator	230 mm	3,5 mm	LAR50-2
Alligator long	230 mm	3,5 mm	LAR50-3
Alligator wide, fenestrated	230 mm	3,5 mm	LAR50-4
Alligator extra wide, fenestrated	230 mm	3,5 mm	LAR50-5



## Cutting forceps, complete matt finish with suction channel




Shape	Lenght	Jaw Diameter	Art. N°
Round cup, Straight	230 mm	2,0 mm	LAR200-2S
Round cup, Upwards	230 mm	2,0 mm	LAR200-2U
Round cup, Right angled	230 mm	2,0 mm	LAR200-2R
Round cup, Left angled	230 mm	2,0 mm	LAR200-2L
Round cup, Straight	230 mm	4,0 mm	LAR200-4S
Round cup, Upwards	230 mm	4,0 mm	LAR200-4U
Round cup, Right angled	230 mm	4,0 mm	LAR200-4R
Round cup, Left angled	230 mm	4,0 mm	LAR200-4L

## Grasping forceps fine serrated, complete matt finish

Shape	Lenght	Jaw Diameter	Art. N°
Fine serrated, Straight	230 mm	2,0 mm	LAR210-2S
Fine serrated, Right angled	230 mm	2,0 mm	LAR210-2R
Fine serrated, Left angled	230 mm	2,0 mm	LAR210-2L
Fine serrated, Straight	230 mm	4,0 mm	LAR210-4S
Fine serrated, Right angled	230 mm	4,0 mm	LAR210-4R
Fine serrated, Left angled	230 mm	4,0 mm	LAR210-4L

That shaft is conically reinforced from distal to proximal end.

Grasping forceps with Laser-Suction channel



Shape	Lenght	Shaft Diameter	Art. N°
Alligator Short	230 mm	3,0 mm	LAR60-1
Alligator long	230 mm	3,5 mm	LAR60-2
Alligator wide, fenestrated	230 mm	3,5 mm	LAR60-3

Rotatable- Laryngo instruments



Shape with conical Tube Shaft 0 3,0 -2.0mm	Lenght	Diameter	Art. N°
Round jaw, right curved	230 mm	2,0 mm	LAR10R-2
Grasping forceps, double action, curved	230 mm	2,0 mm	LAR20R-2
Grasping forceps, triangular, right	230 mm	2,0 mm	LAR28R-2
Grasping forceps, 45° up curved	230 mm	2,0 mm	LAR29R-2
Scissors right curved	230 mm	2,0 mm	LAR30R-2

That shaft is conically reinforced from distal to proximal end.

## Monopolar HF- Instruments



Shape	Lenght	Shaft Diameter	Art. N°
Grasping forceps	230 mm	2,0 mm	LAR80HF-1
Curved	230 mm	2,0 mm	LAR80HF-2
Wide jaw, fenestrated	230 mm	2,0 mm	LAR80HF-3
Oval cup forceps, double action	230 mm	2,0 mm	LAR80HF-4

Complete with isolated shaft, handle and irrigation port  
Demountable in three parts - shaft, handle and insert.

## Flexible grasping and Biopsy Forceps for Laryngo-Fiberscope

Shape	Lenght	Jaw Diameter	Art. N°
Grasping forceps	600 mm	1,0 mm	LAR220-10
Grasping forceps	600 mm	1,7 mm	LAR220-17
Round jaw	600 mm	1,0 mm	LAR221-10
Round jaw	600 mm	1,7 mm	LAR221-17


## SIALENDOSCOPE

### Grasping and Biopsy Forceps for Sialendoscope

Shape	Lenght	Shaft Diameter	Art. N°
Biopsy forceps	270 mm	0,8 mm	BR1108FS-27
Grasping forceps	270 mm	0,8 mm	BR2108FS-27
Biopsy forceps	300 mm	0,8 mm	BR1108FS-30
Grasping forceps	300 mm	0,8 mm	BR2108FS-30



 RG Medical |  RG Technology GmbH Bergstasse 2 97247 Eisenheim Germany

 +49 1525 8601101 / +49 9386-997 997 3

[www.rgmed.com](http://www.rgmed.com)