

MINI MOBILE 12V SURGERY SET

Digital Insufflator, Monitor, Combo Camera 1080p & Light cable

Medical monitors

- ▶ Screen Size 23.8"
- ▶ Resolution 1920 x 1080
- ▶ Brightness 600cd/m²
- ▶ Color Depth 16.7M colors
- ▶ Viewing angle 78/178
- ▶ Response time (GtoG) 14ms
- ▶ Power consumption Typical: 25W



N° Art.
M-23.8CL1

Combo camera LED 300 12V

- ▶ LED 300
- ▶ 1/3" 3.0 Megapixel CMOS Sensor
- ▶ Color: 0.1Lux@F1.2
- ▶ 1920x1080 / 60Hz Full HD
- ▶ 3G SDI 1080p
- ▶ HDMI 1080p
- ▶ C-Video composite video
- ▶ 5300 K
- ▶ Width: 27 cm / Depth: 26 cm / Height: 12 cm



*Full HD Camera + 300w. LED Lightsource
Records in full HD 1080p 30 frames +
Universal light cable

N° Art.
E-CHDI

Digital Insufflator 12V

- ▶ 12 V
- ▶ Width: 26 cm / Depth: 26 cm / Height: 12 cm
- ▶ The highest flow rate available on the market
- ▶ Visual and audible warning signals
- ▶ Safety Relief Valve
- ▶ Internal Leakage Detection
- ▶ Optional gas heating technology



N° Art.
E-IN40