


THE FUTURE IS NOW

NASOPHARYNGOSCOPE **HD**



Instrument complete: S-32HR



Avoid the risk of cross infection with disposable



Multiple connection modes are not subject to the environment and a single device.



The ultra-fine smooth insertion tube makes the patient more "insensitive"



Reusable Handle
Code: E-HR

Disposable: S-32FD
Reposable: S-32FR

Sterilize container

OPTICAL SYSTEM		INSERTION PORTION	
Field of view	100°	Insertion tube OD	3.2mm
Depth of field	3-150mm	Angulation range	U 130° / D 130°
Direction of view	0°	Working length	330 mm ± 10%
		Extension cable	USB, REMO
		Ma. resolution	1280 x 720 p