

# URETHERO-RENOSCOPE FLEXIBLE



- No Cross Infection**
- No Maintenance**
- No Sterilization**
- No optical fiber**
- Handheld, portable, easy use**



Greater Tip  
Range of Motion<sup>TM</sup>



Passive  
Bending




Two-way  
Deflection



1:1 Torque  
Ratio

## 270° ANGULATION



## 12 INCH MONITOR / RECORDER



### SPECIFICATIONS

Field of view	90°
Direction of vision	0° (forward)
Depth of view field	3-50 mm
Image Resolution	1080 p (400x400 native pixels)
Angulation range	Up 270°   Down 270°
Outer diameter insertion tube	7.5Fr / 9.5 Fr
Working Length	650 mm
Working Channel	3.6 Fr
Irrigation	30 mL/min

\*For medical professionals only



RG Medical USA, Inc. | 630 Pleasant Grove Rd, Mount Juliet, TN 37122



+54 9 11 5509-6530

[www.rgmed.com](http://www.rgmed.com)