




# NEW UNIVERSAL MORCESCOPE

**19/22 - 22/24 - 24/26 Fr.**

Only with the RG Morcescope you can choose between MILEP (19/22Fr) or MIBEP (22/24Fr) and can reduce the risk of surgical side effects. After the Laser MILEP or Bipolar MIBEP Enucleation, the Morcescope uses the same sheath for Morcelation avoiding any additional traumatization of the urethra.



**Morcescope** 19 22 24 Fr  
220 mm working length  
5 mm working channel

**Code: U-3-500-00**  
**CMOS Version: U-3-500-00C**

## LASER 19 22 24 Fr. WORKING ELEMENTS




Max. laser fiber  
incl. insulation Ø 1.0 mm

Standard passive

Fr.	Optic	Ref. N°
19 Fr	2.9 mm	U-2-005-19
22 & 24 Fr	4 mm	U-2-005-22

## DUO MONOPOLAR & BIPOLAR WORKING ELEMENT




Working unit, passive  
Standard passive

Fr.	Optic	Ref. N°
19 Fr	2.9 mm	U-1-300-19
22 & 24 Fr	4 mm	U-1-300-22

## SHEATH 19/22 - 22/24 - 24/26 Fr.

19/22Fr.



- 19 Ch./Fr U-1-530-19
- 1 exterior tube U-1-534-19
- 1 inner tube U-1-532-19
- 1 obturator U-1-580-19

22/24 Fr.



- 22 Ch./Fr U-1-530-22
- 1 exterior tube U-1-534-22
- 1 inner tube U-1-532-22
- 1 obturator U-1-580-22

24/26 Fr.



- 24 Ch./Fr U-1-530-24
- 1 exterior tube U-1-534-24
- 1 inner tube U-1-532-24
- 1 obturator U-1-580-24