


NEW UNIVERSAL MORCESCOPE

19/22 - 22/24 - 24/26 Fr.

Only with the RG Morcescope you can choose between MILEP (19/22 Fr) or MIBEP (22/24Fr) and can reducing the risk of surgical side effects. After the Laser MILEP or Bipolar MIBEP Enucleation, the Morcescope uses the same sheath for Morcelation avoiding any additional traumatization of the urethra.



Morcescope 19 22 24 Fr
220 mm working lenght
5 mm working channel

MORCELATOR BLADE


Karl Storz



Stryker



LASER 19 22 24 Fr. WORKING ELEMENTS



Max. laser fiber
incl. insulation Ø 1.0 mm

SHEATH 19/22 - 22/24 - 24/26 Fr.

19/22Fr.



22/24 Fr.



24/26 Fr.

