


DUO ICG 1080p & 4K CAMERA SYSTEM

- ▶ Excellent Real-time Fluorescence Imaging in 1080p or 4k resolution
- ▶ Excellent visualization of visceral perfusion, exposure of biliary anatomy in fluorescent cholangiography & oncological procedures with precise sentinel lymph node mapping
- ▶ Easy Switch:
 - ICG GREEN overlaid Color mode
 - ICG GREEN overlaid Gray Mode
 - White natural color Mode
 - Dual or triple picture - in- picture mode



■ ICG 4Chip 1080 or 4K camera system ■ ICG Duo white LED 300



4K
ULTRAHD

1080p
FULLHD



Excellent Real-time Fluorescence Imaging

Perfect in 1080p or 4k visualization of visceral perfusion, exposure of biliary anatomy in fluorescent cholangiography & oncological procedures with precise sentinel lymph node mapping .

Powerfull 300W LED DUO ICG WHITE

Light for homogeneous, consistently bright endoscopic image. This allows the user to cover two applications with one system – standard white light endoscopy and ICG/NIR fluorescence imaging. Also ideal for Natural color visualization in white light mode

Specifications

- ▶ 1/3" 3.0 Megapixel CMOS Sensor
- ▶ 1920x1080 / 60Hz Full HD HDMI 1080p
- ▶ Touch screen, easy use
- ▶ Compact German Design
- ▶ Width: 27 cm / Depth: 32 cm / Height: 12 cm
- ▶ Ref. Art: E-1080KC

Specifications

- ▶ ICG LED 300
- ▶ White LED 300
- ▶ Touch screen, easy use
- ▶ Compact German Design
- ▶ Width: 27 cm / Depth: 32 cm / Height: 12 cm
- ▶ Ref. Art: E-300L