




THE FUTURE IS NOW

HD RIGID CMOS ENDOSCOPES



Multiple connection modes are not subject to the environment and a single device.



The ultra-fine smooth insertion tube makes the patient more "insensitive"

SPECIFICATIONS		Lenght	Direction of view	110mm	180mm	310mm
Field of view	100°	Ø 3.2 mm	0°	E32-110-00C	E32-180-00C	E32-310-00C
Depth of field	3-150mm		0°	E40-110-00C	E40-180-00C	E40-310-00C
Direction of view	0° -30°	Ø 4.0 mm	30°	E40-110-30C	E40-180-30C	E40-310-30C
Ma. resolution	1280 x 720 p					



RG Medical USA, Inc. | 630 Pleasant Grove Rd, Mount Juliet, TN 37122



615-269-7256

www.rgmed.com