

CF 52 CYSTOSCOPE FLEXIBLE




1280X720P



A variety of connection methods are suitable for various diagnosis and treatment environments



Computer



Pad

SPECIFICATIONS

Field of view	90°
Direction of view	0° (forward)
Depth of view field	3 - 50mm
Angulation range	Up 130° / Down 130°
Insertion tube Outer Diameter	3.2mm
Working length	300mm±10%
Extension cable	USB / REMO
Max. resolution	1280x720 p
Cystoscope flex Ref n°	U-CF-52

**3.2MM
HIGH-DEFINITION
VIDEO ENDOSCOPE**

**ULTRA - THIN
INSERTION**

FULL OF DETAILS

**ULTRA-LIGHT
HANDLE**

MORE STABLE GRIP

**REDUCE THE FATIGUE
OF THE DOCTOR**

**MAKE THE EXAMINATION
SMOOTHER**