




GYNECOLOGY
Hysteroscopy

Developed and manufactured in Germany

■ MEDICAL MONITORS



- ▷ Screen size 24 "
- ▷ Resolution 1920 x 1080
- ▷ Brightness 600cd / m '
- ▷ Color depth
- ▷ 16.7 million colors
- ▷ Viewing angle 78/178
- ▷ Response time (GtoG) 14ms

Code: M-24CL1

■ CAMERAS

1080p wireless camera



- ▷ Without USB
- ▷ Free software
- ▷ Images / Recording

Code: EOC-23

■ DIGITAL INSUFFLATOR



**5 YEAR
WARRANTY**

- ▷ 40 liter digital insufflator
- ▷ Touch screen
- ▷ Day and night mode

Code: E-IN40

1080p USB Camera + 10W LED



- ▷ With USB
- ▷ Free software
- ▷ Images / Recording
- ▷ LED

**Code: EC-USB1
Code: EC-USB1 L**

■ COMBO CAMERA FULL HD

Full HD + 300w camera. LED light source Records in Full HD 1080p 60 frames



- ▷ LED 300
- ▷ 1/3 "CMOS sensor
- ▷ 3.0 megapixels
- ▷ Color: 0.1Lux@F1.2
- ▷ 1920x1080 / 60Hz Full HD
- ▷ 3G SDI 1080p
- ▷ HDMI 1080p
- ▷ Composite video C-Video
- ▷ 5300 K

Code: E-CHD11

* The set includes camera, universal cable and equipment.

■ DUO MONO-BIPOLAR

Endoscopy Ø 2.9 mm



17/18 Fr.

■ Set complet 17/18 Fr	U1-530-17
Inner tube 17 Fr	U1-532-17
Outer tube 18.5 Fr.	U1-534-18
Obturator 17 Fr.	U1-580-17

19/22 Fr.

☐ Set complet 19/22 Fr	U1-530-19
Inner tube 19 Fr	U1-532-19
Outer tube 22 Fr.	U1-534-22
Obturator 19 Fr.	U1-580-19

Endoscopy Ø 4.0 mm



22/24Fr.

■ Set complet 22/24 Fr	U1-530-22
Inner tube 22 Fr	U1-532-22
Outer tube 24 Fr.	U1-534-24
Obturator 22 Fr.	U1-580-22

24/26 Fr.

■ Set complet 24/26 Fr	U1-530-24
Inner tube 24 Fr	U1-532-24
Outer tube 26 Fr.	U1-534-26
Obturator 24 Fr.	U1-580-24

■ ELECTRODES

Cut loop, retrograde 30°



■ 17 Fr.	U1-200-17
□ 19 Fr.	U1-200-19
■ 22 Fr.	U1-400-22
■ 24 Fr.	U1-400-24

Straight cut bow



■ 17 Fr.	U1-202-17
□ 19 Fr.	U1-202-19
■ 22 Fr.	U1-402-22
■ 24 Fr.	U1-402-24

Retrograde cut band



■ 17 Fr.	U1-204-17
□ 19 Fr.	U1-204-19
■ 22 Fr.	U1-404-22
■ 24 Fr.	U1-404-24

Knife electrode



■ 17 Fr.	U1-206-17
□ 19 Fr.	U1-206-19
■ 22 Fr.	U1-406-22
■ 24 Fr.	U1-406-24

Ball electrode 3mm roller



■ 17 Fr.	U1-208-17
□ 19 Fr.	U1-208-19
■ 22 Fr.	U1-408-22
■ 24 Fr.	U1-408-24

Ball electrode 5mm roller



■ 17 Fr.	U1-209-17
□ 19 Fr.	U1-209-19
■ 22 Fr.	U1-409-22
■ 24 Fr.	U1-409-24

■ ELECTRODES BIPOLARES

Cutting loop, retrograde 30°



-  22 Fr. **U1-400-22**
-  24 Fr. **U1-400-24**

Cutting loopstraight



-  22 Fr. **U1-402-22**
-  24 Fr. **U1-402-24**

Knife electrode



-  22 Fr. **U1-406-22**
-  24 Fr. **U1-406-24**

Roller ball electrode 3mm



-  22 Fr. **U1-408-22**
-  24 Fr. **U1-408-24**

Mushroom electrode



-  22 Fr. **U1-416-22**
-  24 Fr. **U1-416-24**

Thick electrode



-  22 Fr. **U1-418-22**
-  24 Fr. **U1-418-24**

Cutting loop, retrograde 0°



-  22 Fr. **U1-403-22**
-  24 Fr. **U1-403-24**

Roller barrel electrode 3mm



-  22 Fr. **U1-414-22**
-  24 Fr. **U1-414-24**

Enucleation electrode



-  22 Fr. **U1-419-22**
-  24 Fr. **U1-419-24**

Push Enucleation electrode



-  22 Fr. **U1-420-22**
-  24 Fr. **U1-420-24**

■ AUTOCLAVE, 134° C / 273° F




Vision angle	Ø 2.0 mm	Ø 2.9 mm	Ø 4.0 mm
0°	E020-240-00	E029-303-00	E040-303-00
12°	E020-240-12	E029-303-12	E040-303-12
30°	E020-240-30	E029-303-30	E040-303-30
70°	E020-240-70	E029-303-70	E040-303-70

■ DIAGNOSTIC SHIRT FOR HYSTEROSCOPY



Scopes	Ext. sheath	Diag. Sheath	Obturator
Ø 2.9 mm	E30-022-05	E30-022-06	E30-022-00
Ø 4.0 mm	E30-023-05	E30-023-06	E30-023-00

■ BETTOCHI SURGICAL SHIRT "STYLE"



Scopes	Working channel Ø	Code
Ø 2.9 mm	5 Charr./Fr	E30-018-05
Ø 4.0 mm	7 Charr./Fr.	E30-018-07




■ FORCEPS



		SEMI RIGID 340 mm	FLEXIBLE 400 mm
Scissor 	5 Ch./Fr 7 Ch./Fr	G30-055-05 G30-055-07	G31-055-05 G30-055-07F
Straight scissors 	5 Ch./Fr 7 Ch./Fr	G30-056-05 G30-056-07	G31-056-05 G31-056-07
Biopsy Clamp 	5 Ch./Fr 7 Ch./Fr	G30-050-05 G30-050-07	G31-050-05 G31-050-07
Aligator Clamp 	5 Ch./Fr 7 Ch./Fr	G30-051-05 G30-051-07	G31-051-05 G31-051-07
Biopsy Grasping 	5 Ch./Fr	G30-052-05	G31-052-05
Tanaculum grasping 	5 Ch./Fr	G30-053-05	G31-053-05
Needle electrode 	5 Ch./Fr	G30-054-05	G31-054-05
Coagulating ball electrode 	5 Ch./Fr	G30-057-05	G31-057-05

■ ELECTRODES



Loop electrode	Code	
5 mm x 5 mm	M90-001	
10 mm x 5 mm	M90-002	
10 mm x 7 mm	M90-003	
10 mm x 10 mm	M90-004	
15 mm x 5 mm	M90-011	
15 mm x 7 mm	M90-012	
15 mm x 8 mm	M90-013	
15 mm x 15 mm	M90-014	
20 mm x 8 mm	M90-021	
20 mm x 10 mm	M90-022	
20 mm x 13 mm	M90-023	
20 mm x 15 mm	M90-024	
20 mm x 20 mm	M90-025	
20 mm x 20 mm	M90-026	
25 mm x 10 mm	M90-027	
Square Electrode		
10 mm x 4 mm	M90-031	
10 mm x 8 mm	M90-032	
Ball Electrode		
3 mm	M90-051	
5 mm	M90-052	

■ LIGHT CABLE



Size	Code
Length 180 cm	E1-035-18

■ VIDEO COLPOSCOPE



Video Colposcope

M99-001

■ TISHER



Tischer
Baby Tischer

G30-070-01
G30-070-02

Kevorkian

G30-070-03

■ MINI 100W LED LIGHT SOURCE



**+ 60%
Light**

■ 100W LED

■ GERMAN DESIGN QUALITY

■ 20.000 + HOURS

■ AUTO-ADJUSTING LIGHT INTENSITY

SPECIFICATIONS

Light Type	100 Watt LED
Ports	Single (Storz), other adapters optional.
Light Intensity Control	Auto Adjust
Color Temp.	6000-6500°K
LED Array Life	20.000 + Hours
Dimensions	9x12x7cm
Weight	280 grs
Warranty	2 Years
Code	E-L100

Limited Warranty: RG Medical instruments are warranted to be free from defects in material and workmanship for the life of the instrument. Any instrument found to be defective will be repaired or replaced at RG Medical's discretion. This warranty is void if the instrument is not used for the surgical purpose for which it was designed or if it is serviced by someone other than RG Medical.

* Please note that the manufacturer only provides a stock image for this product group. The appearance of the product may be different. Specific product details are listed in the description.

■ PORTABLE 10W LED LIGHT



N° Ref: E-L10P

FEATURES

Compact and portable, the built-in lithium battery allows the mini portable cold light source to work for more than 2 hours after recharging.

We design and manufacture it for use with a rigid or flexible endoscope.

A particularly well-conceived safety feature is the automatic switch from the main light source to the spare light source in the event of a light source failure.

Therefore, the light source is always ready for emergency use and contributes to patient safety.



the optician.

Way to use it in

"MINI" LED 10w

N° Ref: E-HL10MP



Includes:

- ▷ LED device
- ▷ 2 rechargeable attitudes
- ▷ Charger
- ▷ Optional small travel case

N° Ref: E-L10P



■ HEADLIGHT SPECIAL

10W ULTRA-BRIGHT & LITE Weight Headlight



*also with focus option

Focusable **E-HL10**
E-HL10F

10W ULTRA-BRIGHT Weight Headlight



10W ULTRA-BRIGHT **E-HL10U**

LED 5W HEADLIGHT Wireless



LED 5W HEADLIGHT Wireless **E-HL5I**

Characteristics

Applications:
Surgical, for dentistry,
veterinary

Lighting type:
LED

Other features:
With rechargeable battery, for binocular
loupe, with portable light source

Headlights include:



- ▶ Portable travel case
- ▶ 1 charger
- ▶ 2 rechargeable lithium-ion batteries



MEDICAL

Technology Germany

 RG Medical USA, Inc. |  630 Pleasant Grove Rd, Mount Juliet, TN 37122
 +1 615 692 1375

www.rgmed.com