





NEPROSCOPIO

5 MM CON DOS CANALES DE FLUIDO CONTINUO DE 4,5MM



Extremo proximal	Longitud de trabajo	Diámetro	Código
 Ocular de 90°	200 mm	5 mm	U4-800-20
	250 mm	5 mm	U4-800-25
 Ocular de 45°	200 mm	5 mm	U4-801-20
	250 mm	5 mm	U4-801-25
 Vaina para nefroscopio con válvula fija	■ 200 mm	24 Fr.	U4-824-20
	■ 250 mm	24 Fr.	U4-824-25
	■ 200 mm	26 Fr.	U4-826-20
	■ 250 mm	26 Fr.	U4-826-25
Nephroscope sheath, rotating stopcock	■ 200 mm	24 Fr.	U4-834-20
	■ 250 mm	24 Fr.	U4-834-25
	■ 200 mm	26 Fr.	U4-836-20
	■ 250 mm	26 Fr.	U4-836-25
 Nephroscope obturator	■ 200 mm	24 Fr.	U4-814-20
	■ 250 mm	24 Fr.	U4-814-25
	■ 200 mm	26 Fr.	U4-816-20
	■ 250 mm	26 Fr.	U4-816-25