

LARYNGOSCOPE

HD

- ▶ Crystal clear image
- ▶ Ergonomic design
- ▶ Sapphire glass lenses

**BRIGHTEST LIGHT
IN THE INDUSTRY!**



LARYNGOSCOPE

300 mm / 2.8 mm	SI-300-2.8
300 mm / 3.0 mm	SI-300-3.0
300 mm / 3.4 mm	SI-300-3.4

with working channel

300 mm / 3.8 mm	SI-300-3.8
-----------------	------------

QUALITY MADE IN GERMANY

FLEXIBLE ENDOSCOPES

Laryngoscope 2.8 / 3.2 / 3.4mm

The new German designed fiber-optic Laryngoscope multiple lens optical system gives a crystal clear image. The extremely thin entry tube is highly flexible and therefore precise work can be done on inaccessible regions. Bigger wide angle with exceptional depth of focus. Illumination of the subject under examination is achieved through two light guides. The flexible entry tube moves precisely to the command of an ergonomically designed well-balanced handpiece. Standard connection is sealed from moisture.


Ergonomic handpiece



Fluent and dynamical move



Camera and Bright Light tip



- Crystal clear image
- Ergonomic design
- Sapphire glass lenses