

PORTABLE 10W LED LIGHT



Features

Compact and portable, the built-in lithium battery allows the mini portable cold light source to work for more than 2 hours after recharging. We design and manufacture it for use with a rigid or flexible endoscope.

A particularly well-conceived safety feature is the automatic switch from the main light source to the spare light source in the event of a light source failure.

Therefore, the light source is always ready for emergency use and contributes to patient safety.



Way to use it in the optician.

“MINI” LED 10W

N° Ref: E-HLI0MP



Includes:

- LED device
- 2 rechargeable attitudes
- Charger
- Optional small travel case

10W LED head-light **E-HLI0P**

