

# CAMERA SYSTEMS

## ■ 4K CAMERA



N° Ref: E-C4k1

### Video Output & Recording 4k



N° Ref: EC-4K

### Specifications

- ▷ CCD: 1/1.8 sony cmos.
- ▷ Cable length: 3 meters.
- ▷ Camera head: waterproof: IPX 67
- ▷ Size camera:  $\phi$ 40\*98mm
- ▷ Ccu output: HDMI+DVI+SDI+USB.
- ▷ Programmable camera buttons
- ▷ 7" touch screen with video recording
- ▷ 5" touch screen without video recording
- ▷ Record: MP4,
- ▷ Capture: JPG

## ■ 1080p HD CAMERA



N° Ref:E-CHD1080

### Video Output & Recording HD 1080



N° Ref: EC-HD

### Specifications

- ▷ 1/3 "CMOS sensor
- ▷ 3.0 megapixels
- ▷ Color: 0.1Lux@F1.2
- ▷ 1920x1080 / 60Hz Full HD
- ▷ 3G SDI 1080p
- ▷ HDMI 1080p
- ▷ Composite video C-Video

## ■ CAMERAS

### Stroboscope USB



N° Ref: EC-USB-S

- ▷ With USB
- ▷ Free software
- ▷ Images / Recording

### 1080p FULL HD



N° Ref: EC-USB1

- ▷ With USB 1080 FULL HD
- ▷ Free software
- ▷ Images / Recording