


2026
Product catalog



■ **GYNECOLOGY**
Hysteroscopy

Developed and manufactured in Germany

■ MEDICAL MONITORS



- ▷ Screen size 24 "
- ▷ Resolution 1920 x 1080
- ▷ Brightness 600cd / m '
- ▷ Color depth
- ▷ 16.7 million colors
- ▷ Viewing angle 78/178
- ▷ Response time (GtoG) 14ms

Code: M-24CL1

■ CAMERAS

1080p wireless camera



- ▷ Without USB
- ▷ Free software
- ▷ Images / Recording

Code: EOC-23

■ DIGITAL INSUFFLATOR



5 YEAR WARRANTY

- ▷ 40 liter digital insufflator
- ▷ Touch screen
- ▷ Day and night mode

Code: E-IN40

1080p USB Camera + 10W LED



- ▷ With USB
- ▷ Free software
- ▷ Images / Recording
- ▷ LED

**Code: EC-USB1
Code: EC-USB1 L**

■ COMBO CAMERA FULL HD

Full HD + 300w camera. LED light source Records in Full HD 1080p 60 frames



- ▷ LED 300
- ▷ 1/3 "CMOS sensor
- ▷ 3.0 megapixels
- ▷ Color: 0.1Lux@F1.2
- ▷ 1920x1080 / 60Hz Full HD
- ▷ 3G SDI 1080p
- ▷ HDMI 1080p
- ▷ Composite video C-Video
- ▷ 5300 K

Code: E-CHD11

* The set includes camera, universal cable and equipment.

■ NEW DUO MONO-BIPOLAR 15, 17 y 19 Fr.

Endoscopy Ø 2.9 mm



15 Fr.

■ Set complet 15 Fr	U1-530-15
Inner tube 15 Fr	U1-532-15
Outer tube 15 Fr.	U1-534-15
Obturator 15 Fr.	U1-580-15
Working element	U1-300-15



17/18 Fr.

■ Set complet 17/18 Fr	U1-530-17
Inner tube 17 Fr	U1-532-17
Outer tube 18.5 Fr.	U1-534-17
Obturator 17 Fr.	U1-580-17
Working element	U1-300-17



19/22 Fr.

□ Set complet 19/22 Fr	U1-530-19
Inner tube 19 Fr	U1-532-19
Outer tube 22 Fr.	U1-534-19
Obturator 19 Fr.	U1-580-19
Working element	U1-300-19

■ DUO MONO-BIPOLAR 22 y 24 Fr.

Endoscopy Ø 4 mm



22/24 Fr.


■ Set complet 22/24 Fr	U1-530-22
Inner tube 22 Fr	U1-532-22
Outer tube 24 Fr.	U1-534-22
Obturator 22 Fr.	U1-580-22
Working element	U1-300-22



24/26 Fr.

■ Set complet 24/26 Fr	U1-530-24
Inner tube 24 Fr	U1-532-24
Outer tube 26 Fr.	U1-534-24
Obturator 24 Fr.	U1-580-24
Working element	U1-300-22

■ WORKING ELEMENT



■ Working unit, passive, Bipolar

Standard passive bipolar

	15 Fr.	U1-300-15
	17 Fr.	U1-300-17
	19 Fr.	U1-300-19
	22 - 24 Fr.	U1-300-22

■ Working unit, active, Bipolar

Standard

	15 Fr.	U1-310-15
	17 Fr.	U1-310-17
	19 Fr.	U1-310-19
	22 - 24 Fr.	U1-310-22

■ ELECTRODES

Cut loop, retrograde 30°



	15 Fr.	U1-200-15
	17 Fr.	U1-200-17
	19 Fr.	U1-200-19
	22 Fr.	U1-400-22
	24 Fr.	U1-400-24

Straight cut bow



	15 Fr.	U1-202-15
	17 Fr.	U1-202-17
	19 Fr.	U1-202-19
	22 Fr.	U1-402-22
	24 Fr.	U1-402-24

Retrograde cut band



	17 Fr.	U1-204-17
	19 Fr.	U1-204-19
	22 Fr.	U1-404-22
	24 Fr.	U1-404-24

Knife electrode



	15 Fr.	U1-206-15
	17 Fr.	U1-206-17
	19 Fr.	U1-206-19
	22 Fr.	U1-406-22
	24 Fr.	U1-406-24

Ball electrode 3mm roller



	17 Fr.	U1-208-17
	19 Fr.	U1-208-19
	22 Fr.	U1-408-22
	24 Fr.	U1-408-24

Ball electrode 5mm roller



	17 Fr.	U1-209-17
	19 Fr.	U1-209-19
	22 Fr.	U1-409-22
	24 Fr.	U1-409-24

■ AUTOCLAVE, 134° C / 273° F




Vision angle	Ø 2.0 mm	Ø 2.9 mm	Ø 4.0 mm
0°	E020-240-00	E029-303-00	E040-303-00
12°	E020-240-12	E029-303-12	E040-303-12
30°	E020-240-30	E029-303-30	E040-303-30
70°	E020-240-70	E029-303-70	E040-303-70

■ DIAGNOSTIC SHIRT FOR HYSTEROSCOPY



Scopes	Ext. sheath	Diag. Sheath	Obturador
Ø 2.9 mm	E30-022-05	E30-022-06	E30-022-00
Ø 4.0 mm	E30-023-05	E30-023-06	E30-023-00









■ BETTOCHI SURGICAL SHIRT "STYLE"



Scopes	Working channel Ø	Code
Ø 2.9 mm	5 Charr./Fr	E30-018-05
Ø 4.0 mm	7 Charr./Fr.	E30-018-07

■ FORCEPS




Scissor 	5 Ch./Fr 7 Ch./Fr
Straight scissors 	5 Ch./Fr 7 Ch./Fr
Biopsy Clamp 	5 Ch./Fr 7 Ch./Fr
Aligator Clamp 	5 Ch./Fr 7 Ch./Fr
Biopsy Grasping 	5 Ch./Fr
Tanaculum grasping 	5 Ch./Fr
Needle electrode 	5 Ch./Fr
Coagulating ball electrode 	5 Ch./Fr

SEMI RIGID 340 mm	FLEXIBLE 400 mm
G30-055-05 G30-055-07	G31-055-05 G31-055-07
G30-056-05 G30-056-07	G31-056-05 G31-056-07
G30-050-05 G30-050-07	G31-050-05 G31-050-07
G30-051-05 G30-051-07	G31-051-05 G31-051-07
G30-052-05	G31-052-05
G30-053-05	G31-053-05
G30-054-05	G31-054-05
G30-057-05	G31-057-05

■ ELECTRODES



Loop electrode	Code
5 mm x 5 mm	M90-001
10 mm x 5 mm	M90-002
10 mm x 7 mm	M90-003
10 mm x 10 mm	M90-004
15 mm x 5 mm	M90-011
15 mm x 7 mm	M90-012
15 mm x 8 mm	M90-013
15 mm x 15 mm	M90-014
20 mm x 8 mm	M90-021
20 mm x 10 mm	M90-022
20 mm x 13 mm	M90-023
20 mm x 15 mm	M90-024
20 mm x 20 mm	M90-025
20 mm x 20 mm	M90-026
25 mm x 10 mm	M90-027
Square Electrode	
10 mm x 4 mm	M90-031
10 mm x 8 mm	M90-032
Ball Electrode	
3 mm	M90-051
5 mm	M90-052



■ LIGHT CABLE



Size	Code
Length 180 cm	E1-035-18

■ VIDEO COLPSCOPE



Video Colposcope	M99-001
------------------	---------

■ TISHER



Tischer	G30-070-01
Baby Tischer	G30-070-02

Kevorkian	G30-070-03
-----------	------------

■ MINI 100W LED LIGHT SOURCE



**+ 60%
Light**

- 100W LED
- GERMAN DESIGN QUALITY
- 20.000 + HOURS
- AUTO-ADJUSTING LIGHT INTENSITY

SPECIFICATIONS

Light Type	100 Watt LED
Ports	Single (Storz), other adapters optional.
Light Intensity Control	Auto Adjust
Color Temp.	6000-6500°K
LED Array Life	20.000 + Hours
Dimensions	9x12x7cm
Weight	280 grs
Warranty	2 Years
Code	E-L100

Limited Warranty: RG Medical instruments are warranted to be free from defects in material and workmanship for the life of the instrument. Any instrument found to be defective will be repaired or replaced at RG Medical's discretion. This warranty is void if the instrument is not used for the surgical purpose for which it was designed or if it is serviced by someone other than RG Medical.

* Please note that the manufacturer only provides a stock image for this product group. The appearance of the product may be different. Specific product details are listed in the description.

■ PORTABLE 10W LED LIGHT



N° Ref: E-L10P

FEATURES

Compact and portable, the built-in lithium battery allows the mini portable cold light source to work for more than 2 hours after recharging.

We design and manufacture it for use with a rigid or flexible endoscope.

A particularly well-conceived safety feature is the automatic switch from the main light source to the spare light source in the event of a light source failure.

Therefore, the light source is always ready for emergency use and contributes to patient safety.



the optician.

Way to use it in

“MINI” LED 10w

N° Ref: E-HL10MP



Includes:

- ▷ LED device
- ▷ 2 rechargeable attitudes
- ▷ Charger
- ▷ Optional small travel case

N° Ref: E-L10P



■ HEADLIGHT SPECIAL

10W ULTRA-BRIGHT & LITE Weight Headlight



*also with focus option

Focusable **E-HL10**
E-HL10F

10W ULTRA-BRIGHT Weight Headlight



10W ULTRA-BRIGHT **E-HL10U**

LED 5W HEADLIGHT Wireless



LED 5W HEADLIGHT Wireless **E-HL5I**

Characteristics

Applications:
Surgical, for dentistry,
veterinary

Lighting type:
LED

Other features:
With rechargeable battery, for binocular
loupe, with portable light source

Headlights include:



- ▶ Portable travel case
- ▶ 1 charger
- ▶ 2 rechargeable lithium-ion batteries



MEDICAL

Technology Germany

 RG Medical USA, Inc. |  630 Pleasant Grove Rd, Mount Juliet, TN 37122

 +1 615 692 1375

www.rgmed.com