



■ LARYNGOSCOPY

Developed and manufactured in Germany



03 Handles



04 Round cup forceps with Handle L-001-01



05 Round cup forceps with Handle AA



06 Oval cup forceps with Handle L-001-01

- Oval cup forceps with Handle AA



07 Suction Tube mountable on the handle

- Triangular with Handle L-001-01



08 Triangular with Handle AA

09 Grasping forceps with Handle L-001-01

10 Grasping forceps with Handle AA



- Suction tube mountable on the handle

11 Scissors with Handle L-001-01

12 Scissors with Handle AA



- 13** Round cup forceps Kleinsasser with Suction channel
- 14** Vocal cord spreader, Needle holder, Shaft curved forceps
- 15** Grasping forceps
- 16** Cutting forceps, complete matt finish with suction channel
 - Grasping forceps fine serrated, complete matt finish
- 17** Grasping forceps with Laser-Suction channel
 - Rotatable - Laryngo instruments
- 18** Monopolar HF - instruments
 - Flexible grasping ans forceps dor laryngo-fiberscope
 - Grasping and biopsy forceps for Slendoscope

2026
Product catalog

www.rgmed.com



Code
L-001-01



Code
L-001-02



Code
L-001-03

Round cup forceps with Handle L-001-01


Shape with conical Tube
Shaft 0.3, 0.2, 0mm

	Lenght	Diameter	Code
Round cup, jaw straight	230 mm	1,0 mm	L1-230-01
Round cup, jaw upwards	230 mm	1,0 mm	L1-230-02
Round cup, jaw right angled	230 mm	1,0 mm	L1-230-03
Round cup, jaw left angled	230 mm	1,0 mm	L1-230-04
Round cup, jaw straight	180 mm	1,0 mm	L1-180-01
Round cup, jaw upwards	180 mm	1,0 mm	L1-180-02
Round cup, jaw right angled	180 mm	1,0 mm	L1-180-03
Round cup, jaw left angled	180 mm	1,0 mm	L1-180-04
Round cup, jaw straight	250 mm	2,0 mm	L1-250-02
Round cup, jaw upwards	250 mm	2,0 mm	L1-250-03
Round cup, jaw right angled	250 mm	2,0 mm	L1-250-04
Round cup, jaw left angled	250 mm	2,0 mm	L1-250-05
Round cup, jaw straight	230 mm	3,0 mm	L1-230-33
Round cup, jaw upwards	230 mm	3,0 mm	L1-230-34
Round cup, jaw right angled	230 mm	3,0 mm	L1-230-35
Round cup, jaw left angled	230 mm	3,0 mm	L1-230-36
Round cup, jaw straight	230 mm	4,0 mm	L1-230-40
Round cup, jaw upwards	230 mm	4,0 mm	L1-230-41
Round cup, jaw right angled	230 mm	4,0 mm	L1-230-42
Round cup, jaw left angled	230 mm	4,0 mm	L1-230-43
Round cup, jaw straight	230 mm	5,0 mm	L1-230-50

Round cup forceps with Handle AA




Shape with conical Tube Shaft 0 3,0 -2.0mm

	Lenght	Diameter	Code
Round cup, jaw straight	230 mm	1,0 mm	LA-230-01
Round cup, jaw upwards	230 mm	1,0 mm	LA-230-02
Round cup, jaw right angled	230 mm	1,0 mm	LA-230-03
Round cup, jaw left angled	230 mm	1,0 mm	LA-230-04
Round cup, jaw straight	180 mm	1,0 mm	LA-180-01
Round cup, jaw upwards	180 mm	1,0 mm	LA-180-02
Round cup, jaw right angled	180 mm	1,0 mm	LA-180-03
Round cup, jaw left angled	180 mm	1,0 mm	LA-180-04
Round cup, jaw straight	250 mm	2,0 mm	LA-250-02
Round cup, jaw upwards	250 mm	2,0 mm	LA-250-03
Round cup, jaw right angled	250 mm	2,0 mm	LA-250-04
Round cup, jaw left angled	250 mm	2,0 mm	LA-250-05
Round cup, jaw straight	230 mm	3,0 mm	LA-230-30
Round cup, jaw upwards	230 mm	3,0 mm	LA-230-31
Round cup, jaw right angled	230 mm	3,0 mm	LA-230-32
Round cup, jaw left angled	230 mm	3,0 mm	LA-230-33
Round cup, jaw straight	230 mm	4,0 mm	LA-230-40
Round cup, jaw upwards	230 mm	4,0 mm	LA-230-41
Round cup, jaw right angled	230 mm	4,0 mm	LA-230-42
Round cup, jaw left angled	230 mm	4,0 mm	LA-230-43
Round cup, jaw straight	230 mm	5,0 mm	LA-230-50

With conical shaft and irrigation port

Oval cup forceps with Handle L-001-01




<u>Shape with conical Tube</u> <u>Shaft 0 3,0 -2.0mm</u>	<u>Lenght</u>	<u>Diameter</u>	<u>Code</u>
Oval cup, jaw straight	230 mm	2,0 mm	LA-230-20
Oval cup, jaw straight	230 mm	3,0 mm	LA-231-30
Oval cup, jaw straight	230 mm	4,0 mm	LA-231-40
Oval cup, jaw straight angled	230 mm	2,0 mm	LA-23-20
Oval cup, jaw straight angled	230 mm	3,0 mm	LA-23-30
Oval cup, jaw straight angled	230 mm	4,0 mm	LA-23-40

Oval cup forceps with Handle AA



<u>Shape with conical Tube</u> <u>Shaft 0 3,0 -2.0mm</u>	<u>Lenght</u>	<u>Diameter</u>	<u>Art. N°</u>
Oval cup, jaw straight	230 mm	2,0 mm	LA-232-20
Oval cup, jaw straight	230 mm	3,0 mm	LA-232-30
Oval cup, jaw straight	230 mm	4,0 mm	LA-232-40
Oval cup, jaw straight angled	230 mm	2,0 mm	LA-233-20
Oval cup, jaw straight angled	230 mm	3,0 mm	LA-233-30
Oval cup, jaw straight angled	230 mm	4,0 mm	LA-233-40


Suction Tube mountable on the handle



The suction tube can be ordered separately

Code
L-001-04

Triangular with Handle L-001-01



Shape with conical Tube
Shaft Ø 3,0 - 2,0mm

Lenght

Diameter

Code

Triangular, jaw straight	230 mm	2,0 mm	LA-234-20
Triangular, jaw right angled	230 mm	2,0 mm	LA-234-21
Triangular, jaw left angled	230 mm	2,0 mm	LA-234-22
Triangular, jaw upwards	230 mm	2,0 mm	LA-234-23
Triangular, jaw straight	230 mm	3,0 mm	LA-234-30
Triangular, jaw right angled	230 mm	3,0 mm	LA-234-31
Triangular, jaw left angled	230 mm	3,0 mm	LA-234-32
Triangular, jaw upwards	230 mm	3,0 mm	LA-234-33


Triangular with Handle AA



Shape with conical Tube Shaft 0 3,0 -2.0mm	Lenght	Diameter	Code
Triangular, jaw straight	230 mm	2,0 mm	LA-235-20
Triangular, jaw right andled	230 mm	2,0 mm	LA-235-21
Triangular, jaw left angled	230 mm	2,0 mm	LA-235-22
Triangular, jaw upwards	230 mm	2,0 mm	LA-235-23
Triangular, jaw straight	230 mm	3,0 mm	LA-236-20
Triangular, jaw right angled	230 mm	3,0 mm	LA-236-21
Triangular, jaw left angled	230 mm	3,0 mm	LA-236-22
Triangular, jaw upwards	230 mm	3,0 mm	LA-236-23

Other Tube diameter upon request available

Grasping forceps with Handle L-001-01



<u>Shape with conical Tube</u>	<u>Lenght</u>	<u>Diameter</u>	<u>Code</u>
<u>Shaft 0 3,0 -2.0mm</u>			
Straight, fine serrated	230 mm	1,0 mm	LA-237-01
Left curved, fine serrated	230 mm	1,0 mm	LA-237-02
Right curved, fine serrated	230 mm	1,0 mm	LA-237-03
Straight, fine serrated	180 mm	1,0 mm	LA-180-01
Left curved, fine serrated	180 mm	1,0 mm	LA-180-02
Right curved, fine serrated	180 mm	1,0 mm	LA-180-03
Straight jaw, fine serrated	230 mm	2,0 mm	LA-238-02
Upwards curved jaw , fine serrated	230 mm	2,0 mm	LA-238-03
Left curved, fine serrated	230 mm	2,0 mm	LA-238-04
Right curved, fine serrated	230 mm	2,0 mm	LA-181-01
Straight, fine serrated	180 mm	2,0 mm	LA-181-02
Left curved, fine serrated	180 mm	2,0 mm	LA-181-03
Right curved, fine serrated	180 mm	2,0 mm	LA-181-04
Straight jaw, fine serrated	230 mm	3,0 mm	LA-239-03
Upwards curved jaw , fine serrated	230 mm	3,0 mm	LA-239-04
Left curved, fine serrated	230 mm	3,0 mm	LA-239-05
Right curved, fine serrated	230 mm	3,0 mm	LA-239-06
Straight, alligator teeth	230 mm	2,0 mm	LA-240-02
Straight, alligator teeth	230 mm	3,0 mm	LA-240-03

Other Tube diameter upon request available

Grasping forceps with Handle AA



<u>Shape with conical Tube Shaft 0 3,0 -2.0mm</u>	<u>Lenght</u>	<u>Diameter</u>	<u>Code</u>
Straight, fine serrated	230 mm	1,0 mm	LA-2341-01
Left curved, fine serrated	230 mm	1,0 mm	LA-2341-02
Right curved, fine serrated	230 mm	1,0 mm	LA-2341-03
Straight, fine serrated	180 mm	1,0 mm	LA-182-01
Left curved, fine serrated	180 mm	1,0 mm	LA-182-02
Right curved, fine serrated	180 mm	1,0 mm	LA-182-03
Straight jaw, fine serrated	230 mm	2,0 mm	LA-242-02
Upwards curved jaw , fine serrated	230 mm	2,0 mm	LA-242-03
Left curved, fine serrated	230 mm	2,0 mm	LA-242-04
Right curved, fine serrated	230 mm	2,0 mm	LA-183-01
Straight, fine serrated	180 mm	2,0 mm	LA-183-02
Left curved, fine serrated	180 mm	2,0 mm	LA-183-03
Right curved, fine serrated	180 mm	2,0 mm	LA-183-04
Straight jaw, fine serrated	230 mm	3,0 mm	LA-243-03
Upwards curved jaw , fine serrated	230 mm	3,0 mm	LA-243-04
Left curved, fine serrated	230 mm	3,0 mm	LA-243-05
Right curved, fine serrated	230 mm	3,0 mm	LA-243-06
Straight, alligator teeth	230 mm	2,0 mm	LA-244-02
Straight, alligator teeth	230 mm	3,0 mm	LA-244-03


Suction Tube mountable on the handle



The suction tube can be ordered separately

Code
L-001-05

Scissors with Handle L-001-01



Shape with conical Tube **Shaft 0 3,0 -2.0mm**

	Lenght	Diameter	Code
Straight	230 mm	1,0 mm	LA-245-01
Left curved	230 mm	1,0 mm	LA-245-02
Right curved	230 mm	1,0 mm	LA-245-03
Straight,	180 mm	1,0 mm	LA-184-01
Left curved	180 mm	1,0 mm	LA-184-02
Right curved	180 mm	1,0 mm	LA-184-03
Straight	250 mm	2,0 mm	LA-251-01
Left curved	250 mm	2,0 mm	LA-251-02
Right curved	250 mm	2,0 mm	LA-251-03
Straight	230 mm	3,0 mm	LA-246-01
Left curved	230 mm	3,0 mm	LA-246-02
Right curved	230 mm	3,0 mm	LA-246-03



Shape	Lenght	Diameter	Code
45° angled curved	230 mm	3,0 mm	LA-045-03

Scissors with Handle AA



Shape with conical Tube Shaft 0 3,0 -2,0mm

	<u>Lenght</u>	<u>Diameter</u>	<u>Code</u>
Straight	230 mm	1,0 mm	AA30-1S
Left curved	230 mm	1,0 mm	AA30-1L
Right curved	230 mm	1,0 mm	AA30-1R
Straight,	180 mm	1,0 mm	AA35-1S
Left curved	180 mm	1,0 mm	AA35-1L
Right curved	180 mm	1,0 mm	AA35-1R
Straight	250 mm	2,0 mm	AA30-2S
Left curved	250 mm	2,0 mm	AA20-2L
Right curved	250 mm	2,0 mm	AA20-2R
Straight	230 mm	3,0 mm	AA30-3S
Left curved	230 mm	3,0 mm	AA30-3L
Right curved	230 mm	3,0 mm	AA30-3R



<u>Shape</u>	<u>Lenght</u>	<u>Diameter</u>	<u>Code</u>
45° angled curved	230 mm	3,0 mm	A45-3U

Round cup forceps Kleinsasser with Suction channel



Shape with conical Tube **Shaft 0 3,0 -2.0mm**

	<u>Lenght</u>	<u>Diameter</u>	<u>Code</u>
Round cup, Straight	230 mm	2,0 mm	LA-247-01
Round cup, Upwards	230 mm	2,0 mm	LA-247-02
Round cup, Left	230 mm	2,0 mm	LA-247-03
Round cup, Right	230 mm	2,0 mm	LA-247-04
Round cup, Straight	230 mm	4,0 mm	LA-248-01
Round cup, Upwards	230 mm	4,0 mm	LA-248-02
Round cup, Left	230 mm	4,0 mm	LA-248-03
Round cup, Right	230 mm	4,0 mm	LA-248-04

Vocal cord spreader, Needle holder, Shaft curved forceps



Code

LA-095-20

Shape with conical Tube
Shaft 0 3,0 -2,0mm

Lenght

Shaft
Diameter

Code

Vocal cord spreader diameter 4,5

250 mm

2,0 mm

LA-070-20

Larynx-, Alligator jaw
shaft 80° downwards

260 mm

2,0 mm

LA-070-21

Larynx-forceps, round jaw 45°
shaft 80° downwards

260 mm

2,0 mm

LA-070-22

Larynx-, Alligator jaw
shaft 80° downwards

260 mm

2,0 mm

LA-070-23

Needle holder forceps with ratchet
handle

250 mm

2,0 mm

LA-070-24



Code

LA-070-20

Grasping forceps



Code
LA-050-01



Code
LA-050-02



Code
LA-050-03



Code
LA-050-04



Code
LA-050-05

Shape with conical Tube **Shaft 0 3,0 -2.0mm**

	Lenght	Diameter	Code
Fine serrated, small	230 mm	2,0 mm	LA-050-01
Alligator	230 mm	3,5 mm	LA-050-02
Alligator long	230 mm	3,5 mm	LA-050-03
Alligator wide, fenestrated	230 mm	3,5 mm	LA-050-04
Alligator extra wide, fenestrated	230 mm	3,5 mm	LA-050-05

Cutting forceps, complete matt finish with suction channel




<u>Shape</u>	<u>Lenght</u>	<u>Jaw Diameter</u>	<u>Code</u>
Round cup, Straight	230 mm	2,0 mm	LA-200-02
Round cup, Upwards	230 mm	2,0 mm	LA-200-03
Round cup, Right angled	230 mm	2,0 mm	LA-200-04
Round cup, Left angled	230 mm	2,0 mm	LA-200-05
Round cup, Straight	230 mm	4,0 mm	LA-201-01
Round cup, Upwards	230 mm	4,0 mm	LA-201-02
Round cup, Right angled	230 mm	4,0 mm	LA-201-03
Round cup, Left angled	230 mm	4,0 mm	LA-201-04

Grasping forceps fine serrated, complete matt finish

<u>Shape</u>	<u>Lenght</u>	<u>Jaw Diameter</u>	<u>Code</u>
Fine serrated, Straight	230 mm	2,0 mm	LA-210-02
Fine serrated, Right angled	230 mm	2,0 mm	LA-210-03
Fine serrated, Left angled	230 mm	2,0 mm	LA-210-04
Fine serrated, Straight	230 mm	4,0 mm	LA-211-02
Fine serrated, Right angled	230 mm	4,0 mm	LA-211-03
Fine serrated, Left angled	230 mm	4,0 mm	LA-211-04

That shaft is conically reinforced from distal to proximal end.

Grasping forceps with Laser-Suction channel



Shape	Lenght	Shaft Diameter	Code
Alligator Short	230 mm	3,0 mm	LA-060-01
Alligator long	230 mm	3,5 mm	LA-060-02
Alligator wide, fenestrated	230 mm	3,5 mm	LA-060-03

Rotatable- Laryngo instruments



<u>Shape with conical Tube</u> <u>Shaft 0 3,0 -2,0mm</u>	Lenght	Diameter	Code
Round jaw, right curved	230 mm	2,0 mm	LA-010-01
Grasping forceps, double action, curved	230 mm	2,0 mm	LA-010-02
Grasping forceps, triangular, right	230 mm	2,0 mm	LA-010-03
Grasping forceps, 45° up curved	230 mm	2,0 mm	LA-010-04
Scissors right curved	230 mm	2,0 mm	LA-010-05

That shaft is conically reinforced from distal to proximal end.

Monopolar HF- Instruments



Shape	Lenght	Shaft Diameter	Code
Grasping forceps	230 mm	2,0 mm	LA-080-01
Curved	230 mm	2,0 mm	LA-080-02
Wide jaw, fenestrated	230 mm	2,0 mm	LA-080-03
Oval cup forceps, double action	230 mm	2,0 mm	LA-080-04

Complete with isolated shaft, handle and irrigation port Demountable in three parts - shaft, handle and insert.

Flexible grasping and Biopsy Forceps for Laryngo-Fiberscope

Shape	Lenght	Jaw Diameter	Code
Grasping forceps	600 mm	1,0 mm	LA-600-10
Grasping forceps	600 mm	1,7 mm	LA-600-17
Round jaw	600 mm	1,0 mm	LA-601-10
Round jaw	600 mm	1,7 mm	LA-601-17

SIALENDOSCOPE

Grasping and Biopsy Forceps for Sialendoscope

Shape	Lenght	Shaft Diameter	Code
Biopsy forceps	270 mm	0,8 mm	LA-270-08
Grasping forceps	270 mm	0,8 mm	LA-270-08
Biopsy forceps	300 mm	0,8 mm	LA-300-08
Grasping forceps	300 mm	0,8 mm	LA-301-08

