



SINOSCOPIA

Desarrollado y fabricado en Alemania.



03

Opticas - compatibilidad

04

Autoclave

05

Morce/Nefroscopio Universal

06

AAO especial

07

ENTsepecial

08

Laringoestroboscopia

09

Mangos

10

Elevadores de tabique

11

Pinza cortante nasal

14

Pinzas circulares

15

Pinza sacabocados antral

16

Pinzas sacabocados para el hueso esfenoides

17


Pinza sacabocados.

2026

Catálogo de productos


www.rgmed.com

ENDOSCOPIOS



Además de nuestros endoscopios, producimos endoscopios comparables a las siguientes marcas bajo pedido

COMPATIBILIDAD





Conector de guía de luz compatible

Las ventajas de los endoscopios RG:

- El sistema óptico es menos sensible a la flexión.
- Mayor flexibilidad en comparación con los endoscopios rígidos estándar
- Mayor robustez (el 70% de descuento en todas las reparaciones se debe a lentes de varilla rotas causadas por tubos doblados)
- Cristal de zafiro a prueba de rayones en el extremo distal y el lado ocular
- Reproducción de alta resolución y contraste
- 1 año de espera contra lentes de varilla rotas

AUTOCLAVE 134° C / 273° F



Sinuscope, Ø 2,7 mm, WL 105mm

Direccion de vista	No-autoclavable	Autoclavable
0°	S60-001-00	S60-001-00A
30°	S60-001-30	S60-001-30A
70°	S60-001-70	S60-001-70A

Sinuscope, Ø 2,7 mm, WL 140mm

Direccion de vista	No - autoclavable	Autoclavable
0°	S60-002-00	S60-002-00A
30°	S60-002-30	S60-002-30A
70°	S60-002-70	S60-002-70A

Sinuscope, Ø 2,7 mm, WL 185mm

Direccion de vista	No - autoclavable	Autoclavable
0°	S60-003-00	S60-003-00A
30°	S60-003-30	S60-003-30A
70°	S60-003-70	S60-003-70A

Sinuscope, Ø 4 mm, WL 185mm

Direccion de vista	No - autoclavable	Autoclavable
0°	S60-003-00	S60-003-00A
30°	S60-003-30	S60-003-30A
70°	S60-003-70	S60-003-70A

AUTOCLAVE 134° C / 273° F



Asas para sinusopes



Mini Sinuscopio, Ø 1,7 mm, Largo 107mm

Direccion de vista	No - autoclavable	Autoclavable
0°	S60-000-00	S60-000-00A

Para endoscopio 0° - 70°

Ø	largo	Art. N°
Ø 2.7 mm	111 mm	S60-000-01
Ø 4.0 mm	110 mm	S60-000-02

MORCE / NEFROSCOPIO UNIVERSAL

19/22 - 22/24 - 24/26 Fr.

Solo con el cistoscopio RG puede elegir entre MILEP (19/22 Fr) o MIBEP (22/24 Fr) y reducir el riesgo de efectos secundarios quirúrgicos. Tras la enucleación láser MILEP o bipolar MIBEP, el cistoscopio utiliza la misma vaina para la morcelación, evitando así cualquier traumatismo adicional en la uretra.




ARTROSCOPIO - SINUSCOPIO - ORL

Dirección de visión	Ø 1.8 mm	Ø 4.0 mm
0°	E29-180-00	E40-180-00
12°	E29-180-12	E40-180-12
30°	E29-180-30	E40-180-30
70°	E29-180-70	E40-180-70



BRONCOSCOPIO - CISTOSCOPIO - HISTEROSCOPIO

Dirección de visión	Ø 3.2 mm	Ø 4.0 mm
0°	E029-303-00	E040-303-00
12°	E029-303-12	E040-303-12
30°	E029-303-30	E040-303-30
70°	E029-303-70	E040-303-70



ESPECIAL AAO

MONITORES MEDICOS



- ▷ Screen Size 23.8"
- ▷ Resolution 1920 x 1080
- ▷ Brightness 600cd/m'
- ▷ Color Depth 16.7M colors
- ▷ Viewing angle 78/178
- ▷ Response time (GtoG) 14ms

M-23.8CLI

COMBO

Full HD Camera + 300w. LED Lightsource
Records in full HD 1080p 60 frames



- ▷ LED 300
- ▷ 1/3" 3.0 Megapixel CMOS Sensor
- ▷ Color: 0.1Lux@F1.2
- ▷ 1920x1080 / 60Hz Full HD
- ▷ 3G SDI 1080p
- ▷ HDMI 1080p
- ▷ C-Video composite video
- ▷ 5300 K
- ▷ Width: 27 cm / Depth: 32 cm / Height: 12 cm


Code: E-CHD1

VIDEO NASOFARINGOSCOPIO

DISEÑO DE MANGO
REUTILIZABLE

REUSABLE SCOPE

CONTROL DE VIDEOS/FOTOS



Instrumento completo: S-32HR

UNIVERSAL LIGHT CABLE



Ø	Working length	Code
3.5mm	1800 mm	E1-035-18
3.5mm	2300 mm	E1-035-23
3.5mm	3000 mm	E1-035-30
3.5mm	1800 mm	E1-048-18
3.5mm	2300 mm	E1-048-23
3.5mm	3000 mm	E1-048-30

With silicone-rubber sheath, autoclavable, waterproof,
very flexible, heat resistant from approx -60° to +180° C

MINI 100W LED LIGTH SOURCE



**+60%
Light**

- ▷ Dimensions: 9x12x7cm
 - ▷ Pound: 280
- Code: E-L100

ESPECIAL ENT

MINI 1080p CAMERA



N° Ref: E-CHD12

Specifications

- ▷ 1/3 "CMOS sensor
- ▷ 3.0 megapixels
- ▷ Color: 0.1Lux@F1.2
- ▷ 1920x1080 / 60Hz Full HD
- ▷ 3G SDI 1080p
- ▷ HDMI 1080p
- ▷ Composite video C-Video

CAMERAS

Stroboscope USB



N° Ref: EC-USB-S

1080p FULL HD



N° Ref: EC-USB1

- ▷ With USB 1080 FULL HD
- ▷ Free software
- ▷ Images / Recording

10W LED HEADLIGHT PORTABLE



*also with zoom option

The headlight package includes:

- ▷ Portable travel suitcase
- ▷ 2 lithium-ion batteries each
Rechargeable batteries
- ▷ Charger

Code: Non- Focousable E-HL10
Focousable E-HL10F



10W LED LIGHT PORTABLE



Thanks to the compact and portable built-in lithium battery, the mini handheld cold light source can work longer 2 hours after charging.

"MINI" LED 10w



Includes:

- ▷ LED device
- ▷ 2 rechargeable batteries
- ▷ Charger
- ▷ Small travel case optional



LARYNGOSTROBOSKOP



FREE!



- 💡 Optimal examination light
- 📋 Durable & Noiseless
- 📷 Stung Sharp
- 🔌 Flexibility and compatibility
- 👍 Pleasant benefits

STROBOSCOPE USB

Technical specifications

Frequency range	70-1200 Hz
Slow motion frequency	fix 0.7 Hz and 1.5 Hz
Microphone control	Airborne / structure-borne sound microphone
Filter oscillations	0 to max. 20ms
Sound level display dB	Luminous display
Light output	Digit height 15mm
Data output	Power LED approx. 4000k
Mains connection	Live sound, sound pressure dB
Power consumption	100 to 240 V - 50/60 Hz
Dimension	W 320, H 100, D 230mm
Weight	2.9 kg

TROCAR



170mm - 6 3/4"

3.5 mm Ø
S60-000-10

5 mm Ø
S60-000-15

5 mm Ø
S60-000-20

ASCENSORES DE TABIQUE



200mm - 8"
S60-101-20



1:2



210mm - 8 1/4"
S60-101-21

1:1

1:2

1:1

CATÉTERES DE EUSTAQUIO



1:2

160mm - 6"


○ 1:1

● 1:1

2.5 mm Ø
S60-100-16

3.0 mm Ø
S60-200-16

ASCENSORES DE TABIQUE



PINZAS DE CORTE NASAL



Schmeden
x= 120mm - 4 3/4"



S60-103-11



S60-103-12



S60-103-13



Tomati
x= 120mm - 4 3/4" S60-104-12




Struycken
x= 120mm - 4 3/4" S60-105-01

0.9 mm


1.2 mm

S60-105-02



Henke
x= 120mm - 4 3/4" S60-106-01

TIJERAS NASALES



Dentado, recto
x= 135 mm S60-107-01




Dentado, derecha
S60-107-01



Dentado, izquierda
S60-107-01

PINZAS DE CORTE NASAL



**Vertical recta,
Acción doble**

y= 3mm
x= 100mm - 4" S60-110-03

**Hor. Recta,
Acción doble**

y= 3mm
x= 100mm - 4" S60-111-03

**Hor. Recta,
acción unica**

y= 3mm
x= 100mm - 4" S60-112-03



**45° vertical,
doble acción**

y= 3mm
x= 100mm - 4" S60-113-03



**45° hor.,
doble acción**

y= 3mm
x= 100mm - 4" S60-114-03

**45° hor.,
acción simple, corta**

y= 3mm
x= 100mm - 4" S60-115-03



**70° vertical,
Acción doble**

y= 3mm
x= 100mm - 4" S60-116-03



**70° hor.,
Acción doble**

y= 3mm
x= 100mm - 4" S60-117-03



**110° hor.,
acción unica**

y= 3mm
x= 100mm - 4" S60-118-03


**70° vertical,
acción simple, corta**

y= 3mm
x= 100mm - 4" S60-119-03

**70° hor.,
acción unica**

y= 3mm
x= 100mm - 4" S60-120-03

PINZAS DE CORTE NASAL



45 ° vertical,
Acción doble

y= 4mm

x= 100mm - 4"

S60-110-04

45 ° vertical,
acción simple, corta corta

y= 4mm

x= 100mm - 4"

S60-112-04



70 ° vertical,
Acción doble

y= 4mm

x= 100mm - 4"

S60-113-04

70 ° vertical,
acción simple, corta corta

y= 4mm

x= 100mm - 4"

S60-115-04

45 ° hor.,
doble acción, corta

y= 4mm

x= 100mm - 4"

S60-111-04

70 ° hor.,
doble acción, corta

y= 4mm

x= 100mm - 4"

S60-114-04



PINZAS CIRCULARES



Derecho
y= 180mm - 7"

Curvado hacia arriba
y= 180mm - 7"

3.5mm Ø


S60-130-03

S60-131-03

4.5mm Ø

S60-130-04

S60-131-04




Con canal ascendente
y= 170mm - 6½"

3.5mm Ø

S60-132-03

Pinza sacabocados antral



x= 100mm
1.5mm x 4.5mm

reves, corte para arriba


S60-140-01

al revés, corte a
la derecha

S60-140-02

al revés, corte a
la izquierda

S60-140-03



HEUWIESER

Mordaza fija curvada 90 °

x= 100mm
S60-141-00

PUNZONES DE HUESO ESFENOIDES



1:1

HAJEK

x= 140 mm - 5½"

4.0mm x 4.0mm

S60-150-00



1:1

HAJEK - KOFLER

x= 140 mm - 5½"

3.5mm x 3.5mm

S60-150-01



1:1

HAJEK

x= 140 mm - 5½"

4.0mm x 4.0mm

S60-151-00




1:1

HAJEK - KOFLER

x= 140 mm - 5½"

3.5mm x 3.5mm

S60-151-01



1:1

HAJEK - KOFLER

El eje gira 360 °

x= 180 mm - 7" / 3.5mm

S60-152-00

HUESO ESFENOIDES PERFORADO



HAJEK - CLAUS

200mm - 8"

S60-160-00



CITELLI

170mm - 6¾"

x=60mm



2.2mm
S60-160-01



3.0mm
S60-160-02



KERRISON

x = 90mm - 3½"
200mm - 8"




4mm x 4mm
S60-170-00



5mm x 5mm
S60-170-01



6mm x 6mm
S60-170-02



FERRIS - SMITH - KERRISON

		x= 150mm - 6"	x= 180mm - 7"
	1mm	S60-170-01	S60-171-01
	2mm	S60-170-02	S60-171-02
	3mm	S60-170-03	S60-171-03
	1mm	S60-172-01	S60-173-01
	2mm	S60-172-02	S60-173-02
	3mm	S60-172-03	S60-173-03
	1mm	S60-174-01	S60-175-01
	2mm	S60-174-02	S60-175-02
	3mm	S60-174-03	S60-175-03

FERRIS - SMITH - KERRISON

Reposapiés delgado

	1mm	S60-176-01
	2mm	S60-176-02
	3mm	S60-176-03

