ENT SPECIAL

MINI CAMERA



CHARACTERISTICS

- 1/3" CMOS sensor
- 3.0 megapixels
- Color: 0.1Lux@F1.2
- 1920x1080 / 60Hz Full HD
- 3G SDI 1080p
- HDMI 1080p
- Composite video C-Video

Code: E-CHD12

10W LED HEADLIGHT PORTABLE



*also with zoom option

THE HEADLIGHT PACKAGE INCLUDES:

- Portable travel suitcase
- 2 lithium-ion batteries each Rechargeable batteries
- Charger

Code: Non- Focousable E-HL10
Focusable E-HL10F



USB CAMERA



USB CAMERA STROBOSCOPE

- Mit USB
- Free software
- Pictures / Recording

Code: EC-USB-S

USB CAMERA 1080p FULL HD

- With USB 1080p HD
- Free software
- Pictures / Recording

Code: EC-USB1

10W LED LIGHT PORTABLE



Thanks to the compact and portable built-in lithium battery, the mini handheld cold light source can work longer 2 hours after charging.

INCLUDES:

- LED device
- 2 rechargeable batteries
- Charger
- Small travel case optional

Code: E-HL10P



RG Medical USA, Inc. | 630 Pleasant Grove Rd, Mount Juliet, TN 37122
+54 9 11 5509-6530
www.rgmed.com